



# HINGMED WBP-02A

## 24 Hour Ambulatory Blood Pressure Monitor



# HINGMED WBP-02A

The Hingmed WBP-02A is a wearable, easy to use 24 hour ambulatory blood pressure monitor.


## FEATURES AND BENEFITS:

- Lightweight, compact and motion tolerant.
- Body position recording for enhanced diagnosis
- USB Rechargeable battery
- BP alarm limits available
- Intelligent pulse wave inflation measurement technology
- Comprehensive report
- Validated and recommended by the ESH and BIHS



# SAVE AND PRINT FINAL REPORT TO PDF:

- Report includes statistical summaries for systolic, diastolic, mean, HR, diff, BP load, asleep dip
- View trends, statistics, detail data, graphs, correlation and hourly averages
- Local workstation or server / client installation options
- Save and print to PDF



## INCLUDED IN THE BOX:



## SPECIFICATIONS:

### BLOOD PRESSURE CUFF SIZES

Small Adult cuff (18-26cm)

Middle Adult cuff (22-32cm)

Large Adult Cuff (22-36cm)

Extra Large Adult cuff (33-43cm)

Cotton sleeve liner - one size fits all

### SPECS

Lightweight 105g

USB Rechargeable Li Battery

Low power consumption - one full charge supports 200 measurements

Size 119 x 52 x 21mm

### RANGE

Systolic 40-260mmHg

Diastolic 20-210mmHg

Heart Rate 40-200 Bmp

Accuracy +/- 3mm Hg